

Aurin Arterial Catheter

The AURIN Arterial Catheters are designed for continuous invasive blood pressure monitoring and repeated arterial blood sampling in critical care and surgical environments. Made using high-quality biocompatible materials, these catheters provide reliable performance, smooth insertion and enhanced patient safety.

Continuous BP Monitoring

Repeated Blood Sampling

Enhanced Patient Safety



Features

- High-quality biocompatible materials
- Smooth insertion with enhanced patient safety
- Reliable performance in critical care and surgical settings
- Available as a complete sterile kit

Technical Specifications

Gauge	20G / 22G / 24G
Length	3–5 cm
Kit May Include	Arterial Catheter, Introducer Needle, Guidewire (optional), Injection Cap, Securement Accessories, Sterile Dressing Components

Available SKUs & Ordering Information

7 variants

ITEM CODE	DESCRIPTION	QTY/MASTER CASE	QTY/INNER BOX	SHELF LIFE
90030101	Arterial Catheter Kit 22G/2*4CM (PTFE)	200	20	5 years
90030102	Arterial Catheter Kit 22G/2*6CM (PTFE)	200	20	5 years
90030103	Arterial Catheter Kit 20G/ 3*8CM (PU)	100	10	5 years
90030104	Arterial Catheter Kit 18G/ 4*10CM (PU)	100	10	5 years
90030201	Artylyne Arterial Catheter Kit 16G/5*11CM (PU)	40	5	5 years
90030202	Artylyne Arterial Catheter Kit 16G/5*15CM (PU)	40	5	5 years
90030203	Artylyne Arterial Catheter Kit 16G/5*20CM (PU)	40	5	5 years