


Aurin AV Fistula Needle

The AURIN AVF Needles are specially designed for safe and efficient hemodialysis access. Made using ultra-thin wall stainless steel needles with a siliconized surface for smooth cannulation and optimized blood flow with back eye. An ultra-sharp three-bevel needle tip with back eye enables painless insertion, and colour-coded wings allow quick gauge identification.

Ultra-Sharp Three-Bevel Tip

Siliconized Surface

Colour-Coded Wings



Features

- Ultra-thin wall stainless steel for optimized blood flow
- Colour-coded wings for quick gauge identification
- Three-bevel tip with back eye for painless insertion
- Available in Fixed / Rotatable wings, Single or Twin packing

Technical Specifications

Gauge	15G / 16G / 17G
Needle Length	25 mm
Tube Length	300 mm
Wing Type	Fixed / Rotatable
Packing	Single / Twin

Available SKUs & Ordering Information

6 variants

ITEM CODE	DESCRIPTION	QTY/MASTER CASE	QTY/INNER BOX	SHELF LIFE
90080101	AVF NEEDLE 15G 300MM, FIXED WING , WITH BACKEYE	1000	100	3 years
90080102	AVF NEEDLE 16G 300MM, FIXED WING , WITH BACKEYE	1000	100	3 years
90080103	AVF NEEDLE 17G 300MM, FIXED WING , WITH BACKEYE	1000	100	3 years
90080201	AVF NEEDLE 15G 300MM, ROTATEING WING , WITH BACKEYE	1000	100	3 years
90080202	AVF NEEDLE 16G 300MM, ROTATEING WING , WITH BACKEYE	1000	100	3 years
90080203	AVF NEEDLE 17G 300MM, ROTATEING WING , WITH BACKEYE	1000	100	3 years