


Aurin Central Venous Catheter (CVC)

The AURIN Central Venous Catheter (CVC) kits are designed to provide reliable vascular access for intensive monitoring and complex therapeutic infusions. Made from biocompatible polyurethane, the catheter remains firm during insertion for effortless placement, then softens at body temperature to minimize vascular trauma.

Biocompatible Polyurethane

Softens at Body Temperature

Complete Kit



Features

- Firm during insertion for effortless placement
- Made from biocompatible polyurethane
- Softens at body temperature to minimize vascular trauma
- Supplied as a complete procedural kit

Technical Specifications

Material	Biocompatible Polyurethane
Kit Components	CVC, J-Tip Guidewire, Introducer Needle, Vessel Dilator, Injection Site Cap
Behaviour	Firm on insertion, softens at body temperature

Available SKUs & Ordering Information

7 variants

ITEM CODE	DESCRIPTION	QTY/MASTER CASE	QTY/INNER BOX	SHELF LIFE
90100101	SLCVC 5 FR * 15 /20cm	30	5	5 years
90100102	SLCVC 6 FR * 15/20cm	30	5	5 years
90100103	SLCVC 7FR * 15 /20cm	30	5	5 years
90100111	DLCVC 7 FR * 15cm ANTIMICROBIAL	30	5	5 years
90100112	DLCVC 7 FR * 20cm ANTIMICROBIAL	30	5	5 years
90100121	TLCVC 7 FR * 15cm ANTIMICROBIAL	30	5	5 years
90100122	TLCVC 7 FR * 20cm ANTIMICROBIAL	30	5	5 years