

Aurin Dialyzer

The AURIN dialyzers are engineered with cutting-edge membrane technology to deliver superior clearance of low and middle molecules while uncompromisingly safeguarding patient safety. Designed for maximum biocompatibility and consistent performance, they bring clinical efficiency into perfect alignment.

- High Biocompatibility
- Consistent Clearance
- Low & High Flux



Features

- Advanced Polyethersulfone (PES) / Polysulfone membrane
- Polycarbonate housing with Polyurethane potting
- Available in Low-Flux and High-Flux series
- Gamma-ray sterilized for reliable performance

Technical Specifications

Membrane	Advanced Polyethersulfone (PES) / Polysulfone
Effective Surface Area	Low-Flux 1.3 m ² · High-Flux 1.4 m ²
UF Coefficient (Kuf)	Low-Flux 10–15 · High-Flux 40–65 mL/h/mmHg
Clearance (Urea @ Qb = 200 mL/min)	Low-Flux > 185 · High-Flux > 195 mL/min
Housing & Potting	Polycarbonate / Polyurethane
Sterilization Method	Gamma Ray
Maximum TMP	Low-Flux 500 mmHg · High-Flux 600 mmHg

Available SKUs & Ordering Information

21 variants

ITEM CODE	DESCRIPTION	QTY/MASTER CASE	SHELF LIFE
90090101	DIALYZER PES PC 1.0 SQR MTR	30	3 years
90090102	DIALYZER PES PC 1.3 SQR MTR	30	3 years
90090103	DIALYZER PES PC 1.4 SQR MTR	30	3 years
90090104	DIALYZER PES PC 1.5 SQR MTR	30	3 years
90090105	DIALYZER PES PC 1.6 SQR MTR	30	3 years
90090106	DIALYZER PES PC 1.7 SQR MTR	30	3 years
90090107	DIALYZER PES PC 1.8 SQR MTR	30	3 years
90090108	DIALYZER PES PC 1.9SQR MTR	30	3 years
90090109	DIALYZER PES PC 2.0 SQR MTR	30	3 years
90090201	DIALYZER PES PP 1.0 SQR MTR	30	3 years
90090202	DIALYZER PES PP 1.3 SQR MTR	30	3 years
90090203	DIALYZER PES PP 1.4 SQR MTR	30	3 years
90090204	DIALYZER PES PP 1.5 SQR MTR	30	3 years
90090205	DIALYZER PES PP 1.6 SQR MTR	30	3 years
90090206	DIALYZER PES PP 1.7 SQR MTR	30	3 years
90090207	DIALYZER PES PP 1.8 SQR MTR	30	3 years
90090208	DIALYZER PES PP 1.9SQR MTR	30	3 years
90090209	DIALYZER PES PP 2.0 SQR MTR	30	3 years
90090301	HF DIALYZER PES PC 1.5 SQR MTR	30	3 years
90090302	HF DIALYZER PES PC 1.8 SQR MTR	30	3 years
90090303	HF DIALYZER PES PC 2.0 SQR MTR	30	3 years