


Aurin Hemodialysis Blood Tubing Set

The AURIN Hemodialysis Blood Tubing Set is an essential, single-use sterile extracorporeal circuit that guarantees safe and efficient blood transfer during dialysis treatment. Designed with precision, it features both arterial and venous bloodlines, ensuring uninterrupted, continuous and contamination-free blood flow for optimal patient safety and care.

Zero Leak

Superior Flow Performance

Maximum Safety



Features

- Precision-engineered Luer Lock connectors ensure leak-proof connections
- Medical-grade DEHP-free / non-toxic PVC material, pyrogen-free
- Smooth inner surface reduces clotting & turbulence
- Kink-resistant, flexible tubing suited to modern dialysis centers

Technical Specifications

Material	Medical-grade PVC (DEHP Free)
Sterilization	ETO
Tubing Type	Kink-resistant, flexible
Connector	Leak-proof Luer Lock
Sizes	Adult / Pediatric
Packaging	Individually packed, sterile
Usage	Single-use
Compatibility	Universal
Variant	Post-Pump / Pre-Pump

Available SKUs & Ordering Information

8 variants

ITEM CODE	DESCRIPTION	QTY/MASTER CASE	SHELF LIFE
90070101	BLOOD TUBING SET UNIVERSAL POST PUMP (7 X 10)	30	3 years
90070102	BLOOD TUBING SET UNIVERSAL POST PUMP WITH TP (7 X 10)	30	3 years
90070103	BLOOD TUBING SET UNIVERSAL POST PUMP WITH TP AND VENUS FILTER (7 X 10)	30	3 years
90070104	BLOOD TUBING SET UNIVERSAL POST PUMP WITH TP, VENUS FILTER, DRAINAGE BAG (7 X 10)	30	3 years
90070201	BLOOD TUBING SET UNIVERSAL PRE PUMP (7 X 10)	30	3 years
90070202	BLOOD TUBING SET UNIVERSAL PRE PUMP WITH TP (7 X 10)	30	3 years
90070203	BLOOD TUBING SET UNIVERSAL PRE PUMP WITH TP AND ARTERIAL FILTER (7 X 10)	30	3 years
90070204	BLOOD TUBING SET UNIVERSAL PRE PUMP WITH TP, VENUS FILTER, DRAINAGE BAG (7 X 10)	30	3 years