

Aurin Long Term Hemodialysis Catheter

The AURIN long-term Hemodialysis Catheter is specifically designed to provide reliable vascular access for patients undergoing prolonged hemodialysis treatment. Made from biocompatible medical-grade polyurethane (PU), it offers patient comfort, kink resistance and high flow performance. The butterfly valve with a peel-away design prevents blood leakage during insertion.

Medical-Grade Polyurethane

Step-Tip Flow Design

Peel-Away Butterfly Valve



Features

- Soft medical-grade polyurethane with high-flow lumen design
- Step tip optimizes venous flow dynamics
- Colour-coded clamps, smooth polyester cuffs, rotating suture wing
- Radiopaque body for enhanced visibility
- Kink-resistant tubing with clear extension tube
- Sterile and pyrogen-free

Clinical benefits. Optimal blood flow efficiency for dialysis, reduced thrombosis risk and improved long-term patency enhance patient comfort. The step tip optimizes venous flow dynamics.

Technical Specifications

Material	Medical-grade Polyurethane (PU)
French Sizes	8 Fr to 14.5 Fr
Lengths	18 cm to 60 cm
Tip	Step tip
Valve	Butterfly valve, peel-away design
Sterility	Sterile, pyrogen-free

Available SKUs & Ordering Information

80 variants

ITEM CODE	DESCRIPTION	QTY/MASTER CASE	QTY/INNER BOX	SHELF LIFE
CATHETER				
90010101	DLC 11.5 FR * 13.5 cm	30	5	5 years
90010102	DLC 11.5 FR * 15cm	30	5	5 years
90010103	DLC 11.5 FR * 16cm	30	5	5 years
90010104	DLC11.5 FR * 20cm	30	5	5 years
90010105	DLC11.5 FR * 24cm	30	5	5 years
90010106	DLC 12 FR * 13.5cm	30	5	5 years
90010107	DLC 12 FR * 15 cm	30	5	5 years
90010108	DLC 12 FR * 16cm	30	5	5 years
90010109	DLC 12 FR * 20cm	30	5	5 years
90010110	DLC 12 FR * 24 cm	30	5	5 years
90010201	DLS 11.5 FR * 13.5cm	30	5	5 years
90010202	DLS 11.5 FR * 15 cm	30	5	5 years
90010203	DLS 11.5 FR * 16cm	30	5	5 years
90010204	DLS 11.5 FR * 20 cm	30	5	5 years
90010205	DLS 11.5 FR * 24 cm	30	5	5 years
90010206	DLS 12 FR * 13.5cm	30	5	5 years
90010207	DLS 12 FR * 15cm	30	5	5 years
90010208	DLS 12 FR * 16cm	30	5	5 years
90010209	DLS 12 FR * 20cm	30	5	5 years
90010210	DLS 12 FR * 24cm	30	5	5 years
90020101	TLC 11.5 FR * 13.5 cm	30	5	5 years
90020102	TLC 11.5 FR * 15cm	30	5	5 years
90020103	TLC 11.5 FR * 16cm	30	5	5 years
90020104	TLC11.5 FR * 20cm	30	5	5 years
90020105	TLC11.5 FR * 24cm	30	5	5 years
90020106	TLC 12 FR * 13.5cm	30	5	5 years
90020107	TLC 12 FR * 15 cm	30	5	5 years
90020108	TLC 12 FR * 16cm	30	5	5 years
90020109	TLC 12 FR * 20cm	30	5	5 years
90020110	TLC 12 FR * 24 cm	30	5	5 years
90020201	TLS 11.5 FR * 13.5cm	30	5	5 years
90020202	TLS 11.5 FR * 15 cm	30	5	5 years
90020203	TLS 11.5 FR * 16cm	30	5	5 years

ITEM CODE	DESCRIPTION	QTY/MASTER CASE	QTY/INNER BOX	SHELF LIFE
90020204	TLS 11.5 FR * 20 cm	30	5	5 years
90020205	TLS 11.5 FR * 24 cm	30	5	5 years
90020206	TLS 12 FR * 13.5cm	30	5	5 years
90020207	TLS 12 FR * 15cm	30	5	5 years
90020208	TLS 12 FR * 16cm	30	5	5 years
90020209	TLS 12 FR * 20cm	30	5	5 years
90020210	TLS 12 FR * 24cm	30	5	5 years
COMBO KIT				
90010301	DLC 11.5 FR * 13.5 cm	30	5	5 years
90010302	DLC 11.5 FR * 15cm	30	5	5 years
90010303	DLC 11.5 FR * 16cm	30	5	5 years
90010304	DLC11.5 FR * 20cm	30	5	5 years
90010305	DLC11.5 FR * 24cm	30	5	5 years
90010306	DLC 12 FR * 13.5cm	30	5	5 years
90010307	DLC 12 FR * 15 cm	30	5	5 years
90010308	DLC 12 FR * 16cm	30	5	5 years
90010309	DLC 12 FR * 20cm	30	5	5 years
90010310	DLC 12 FR * 24 cm	30	5	5 years
90010401	DLS 11.5 FR * 13.5cm	30	5	5 years
90010402	DLS 11.5 FR * 15 cm	30	5	5 years
90010403	DLS 11.5 FR * 16cm	30	5	5 years
90010404	DLS 11.5 FR * 20 cm	30	5	5 years
90010405	DLS 11.5 FR * 24 cm	30	5	5 years
90010406	DLS 12 FR * 13.5cm	30	5	5 years
90010407	DLS 12 FR * 15cm	30	5	5 years
90010408	DLS 12 FR * 16cm	30	5	5 years
90010409	DLS 12 FR * 20cm	30	5	5 years
90010410	DLS 12 FR * 24cm	30	5	5 years
90020301	TLC 11.5 FR * 13.5 cm	30	5	5 years
90020302	TLC 11.5 FR * 15cm	30	5	5 years
90020303	TLC 11.5 FR * 16cm	30	5	5 years
90020304	TLC11.5 FR * 20cm	30	5	5 years
90020305	TLC11.5 FR * 24cm	40	5	5 years
90020306	TLC 12 FR * 13.5cm	40	5	5 years
90020307	TLC 12 FR * 15 cm	40	5	5 years
90020308	TLC 12 FR * 16cm	30	5	5 years

ITEM CODE	DESCRIPTION	QTY/MASTER CASE	QTY/INNER BOX	SHELF LIFE
90020309	TLC 12 FR * 20cm	30	5	5 years
90020310	TLC 12 FR * 24 cm	30	5	5 years
90020401	TLS 11.5 FR * 13.5cm	30	5	5 years
90020402	TLS 11.5 FR * 15 cm	40	5	5 years
90020403	TLS 11.5 FR * 16cm	40	5	5 years
90020404	TLS 11.5 FR * 20 cm	40	5	5 years
90020405	TLS 11.5 FR * 24 cm	30	5	5 years
90020406	TLS 12 FR * 13.5cm	30	5	5 years
90020407	TLS 12 FR * 15cm	30	5	5 years
90020408	TLS 12 FR * 16cm	30	5	5 years
90020409	TLS 12 FR * 20cm	40	5	5 years
90020410	TLS 12 FR * 24cm	40	5	5 years