

Aurin Peripherally Inserted Central Catheter (PICC)

The AURIN PIC catheters are designed for safe, reliable and long-term venous access for administration of medications, fluids, parenteral nutrition and blood sampling. Made using PU material, these catheters offer superior flexibility, patient comfort and optimal flow performance.

Long-Term Access

Polyurethane Flexibility

Multi-Lumen Options



Features

- Made using PU material for superior flexibility
- Optimal flow performance and patient comfort
- Suited to medications, fluids, parenteral nutrition and blood sampling
- Available in Single, Double and Triple Lumen configurations

Technical Specifications

Single Lumen	3 Fr & 4 Fr
Double Lumen	4 Fr & 5 Fr
Triple Lumen	5.5 Fr
Lengths	40 cm to 60 cm
Material	Polyurethane (PU)

Available SKUs & Ordering Information

13 variants

ITEM CODE	DESCRIPTION	QTY/MASTER CASE	QTY/INNER BOX	SHELF LIFE
SINGLE-LUMEN PICC LINES				
90040101	P-Lyne 3 Fr x 60 cm Standard PICC Line Catheter	5	1	4 years
90040102	P-Lyne 4 Fr x 60 cm Standard PICC Line Catheter	5	1	4 years
90040103	P-Lyne 5 Fr x 60 cm Standard PICC Line Catheter	5	1	4 years
90040104	P-Lyne 4 Fr x 60 cm Power Injector-CT PICC Line Catheter	5	1	4 years
90040105	P-Lyne 5 Fr x 60 cm Power Injector-CT PICC Line Catheter	5	1	4 years
DOUBLE-LUMEN PICC LINES				
90040201	P-Lyne 4 Fr x 60 cm Standard PICC Line Catheter	5	1	4 years
90040202	P-Lyne 5 Fr x 60 cm Standard PICC Line Catheter	5	1	4 years
90040203	P-Lyne 6 Fr x 60 cm Standard PICC Line Catheter	5	1	4 years
90040204	P-Lyne 4 Fr x 60 cm Power Injector-CT PICC Line Catheter	5	1	4 years
90040205	P-Lyne 5 Fr x 60 cm Power Injector-CT PICC Line Catheter	5	1	4 years
90040206	P-Lyne 6 Fr x 60 cm Power Injector-CT PICC Line Catheter	5	1	4 years
TRIPLE-LUMEN PICC LINES				
90040301	P-Lyne 6 Fr x 60 cm Standard PICC Line Catheter	5	1	4 years
90040301	P-Lyne 6 Fr x 60 cm Power Injector-CT PICC Line Catheter	5	1	4 years