

Aurin Transducer Protector

The AURIN Transducer Protector is a sterile, single-use hydrophobic filter designed to prevent infection during dialysis and to monitor hemodialysis tubing pressure. It effectively blocks blood, fluids and microbial contamination while ensuring accurate pressure transmission. With leak-proof Luer Lock connectors and advanced membrane technology, it offers zero leakage and reliable protection.

Zero Leak

Superior Air Barrier

Infection Control Barrier



Features

- Precision sealing prevents fluid or air leakage under pressure
- Acts as a critical infection control barrier
- Universal compatibility, individually packed and sterile
- Hydrophobic membrane allows air passage but blocks liquids
- Transparent housing for visual inspection

Technical Specifications

Filter Type	Hydrophobic membrane
Connector	Leak-proof Luer Lock
Housing	Transparent for visual inspection
Compatibility	Universal
Packaging	Individually packed, sterile
Usage	Single-use

Available SKUs & Ordering Information

1 variant

ITEM CODE	DESCRIPTION	QTY/MASTER CASE	QTY/INNER BOX	SHELF LIFE
90120101	TRANSDUCER PROTECTOR	1000	100	5 years